

Styletmonitor

Model:
ASM-100

Registration Number: 13B2X00217000003

For your casual use, we offer an affordable price.

Wireless System V1.0

Recharge for up to
1,000 times use.

6mm ϕ 5mm ϕ
CCD camera
LED
Tip of Stylet

STYLET
Model: SL-1

Reusable up to
100 times

2.5-inch
Monitor


LCD Monitor
Model: LCD-M
Transmitter Included

Reception Module
Model: WVM-01

Receives images transmitted from
the LCD monitor and displays them
on the TV monitor.

TV Monitor (not included)

A complete set of our StyletMonitor Model: ASM-100




Stylet

Model: SL-1 Dimension: 370 mm (L)



LCD Monitor

Model: LCD-M (transmitter included)
Dimensions: 73 (W) × 57 (H) × 23 (D) mm



Reception Module

Model: WVM-01
Dimensions: 50 (W) × 13 (H) × 60 (D) mm

〈 In the Box 〉

	Quantity	Model
Stylet	1	SL-1
LCD Monitor	1	LCD-M
Reception Module	1	WVM-01
Reception AC Adaptor	1	WVM-AC
Charger for Monitor	1	LCD-C


PRICE

A set of StyletMonitor
ASM-100 **JPY 170,000**

*An alternative set which includes only Stylet and LCD Monitor is also available.

Example of Use:

to aid tracheal intubation,
visual assessment of inner
ear, etc.



Options

9-inch Monitor



9-inch Monitor with Stand



Accessories



Reception AC Adaptor
Model: WVM-AC



Charger for Monitor
Model: LCD-C



STAR MEDICAL, INC.

5F Taiyo-Building 5-14-12 Nishinippori Arakawa-ku Tokyo 116-0013 Japan TEL:+81-3-5604-1641 FAX:+81-3-5604-1645
E-mail:starmedical@starmedical.co.jp <http://www.starmedical.co.jp>